

HEALTH PROTECTION AGENCY
MICROBIOLOGY SERVICES COLINDALE
VIRUS REFERENCE DEPARTMENT
STANDARD OPERATING PROCEDURE

TITLE: INFLUENZA MUNANA NEURAMINIDASE ACTIVITY AND INHIBITION ASSAY (FLUORESCENT IC₅₀ ASSAY)

SOP NO. V-6815/01-10

AUTHOR: ANGIE LACKENBY

This SOP is provided for information only and should not be distributed further.

Please acknowledge the Respiratory Virus Unit, Health Protection Agency when publishing data generated using this SOP.

SUMMARY

This SOP describes the method to determine influenza virus neuraminidase (NA) activity and sensitivity to neuraminidase inhibitors (NI) using an enzyme assay with a fluorescent substrate. NA activity and NI sensitivity can be determined using the fluorogenic substrate, MUNANA (2' 2'-(4-Methylumbelliferyl)- α -D-N-acetylneuraminic acid sodium salt hydrate). This substrate is cleaved by NA to yield free 4-methylumbelliferone, and the quantitative increase in fluorescence gives a measure of NA activity. The concentration of drug needed for inhibition of enzyme activity by 50% (IC₅₀) is determined by assay in the presence of NIs.

SAFETY

This assay is suitable for tissue culture and egg grown influenza A and influenza B viruses

Good Laboratory Practice and refer to VRD Safety Manual (V6764)

COSHH Risk Assessment No.s:

VB37 Influenza Virus Neuraminidase Inhibition

VB551 Influenza Virus Neuraminidase Inhibition with strains to be handled at CL3

VC893 Influenza virus Neuraminidase Inhibition with Fluorescent substrate

VB545 Control measures for work with pandemic influenza H1N1 viruses

SOP NO. V6815/01-10

COPY NO.

PAGE 1 OF 15

1.0 CROSS REFERENCE

- 1.1 VW0771 IC50 Data Management Workflow
- 1.2 VW0770 IC50 Dilution Ranges and Layouts
- 1.3 VW0772 IC50 Reference Virus Validation Limits
- 1.4 VW0574 Munana Neuraminidase Activity Assay Worksheet
- 1.5 VW0575 IC50 Results Record Sheet: Neuraminidase Inhibitor Susceptibility
- 1.6 SOP V6816 IC50 Outlier Identification and Trend Monitoring
- 1.7 SOP V5402 Haemagglutination (HA) Test for Influenza Virus

2.0 PERSONNEL

- 2.1 All medical microbiologists, clinical scientists, biomedical scientists, and healthcare scientists with suitable training.

3.0 EQUIPMENT

- 3.1 Fluorescence plate reader (355nm and 460nm filters) (useful note 9.1)
- 3.2 Plate shaker
- 3.3 Single channel pipettes suitable for 10 to 900µl volumes and filtered tips
- 3.4 8 and/or 12-well multi-channel pipettes suitable for 10 to 150µl volumes and filtered tips
- 3.5 Multi-well reservoirs (Thermo Electron Cat. No. RTP/08200-10)
- 3.6 10 and 25 ml disposable pipettes and pipette boy
- 3.7 Warm Room (+37°C), Fridge (+4°C), Freezer (- 20°C and - 80°C)
- 3.8 Black 96 well flat bottom plates (Corning 3915 or equivalent)
- 3.9 Adhesive plate sealers (useful note 9.2)

4.0 REAGENTS

- 4.1 Influenza virus isolates (tissue culture/egg fluids) HA ≥16 Units (useful note 9.3)
- 4.2 Subtype matched reference viruses (see section 8 for details)
- 4.3 2-Morpholinoethanesulfonic acid (MES) (Sigma-Aldrich M3671 or equivalent)
- 4.4 Calcium chloride (VWR 5701 or equivalent)
- 4.5 Oseltamivir Carboxylate (Roche. Product no. GS4071 or Ro64-0802)
- 4.6 Zanamivir (Glaxo-Smithkline Product no. GR121167X or GG167)
- 4.7 MUNANA (2' 2'-(4-Methylumbelliferyl)-α-D-N-acetylneuraminic acid sodium salt hydrate) (Sigma-Aldrich M8639)
- 4.8 4-Methylumbelliferone sodium salt (Sigma-Aldrich M1508)
- 4.9 Glycine (VWR 1517 or equivalent)
- 4.10 Absolute Ethanol (VWR 101077Y or equivalent)
- 4.11 Sodium Hydroxide (VWR 101182 or equivalent)
- 4.12 Distilled water

5.0 PREPARATION OF BUFFERS AND SOLUTIONS

All solutions and buffers should be stored at room temperature unless otherwise stated. Working solutions for use in the assay are prepared from master stock solutions where stated, for accuracy. The working solution of MUNANA must be made freshly for each assay.

5.1 Master Stock Solutions and Buffers

| | |
|--------------------------------------|---|
| 325mM MES: | 31.72g MES in 500ml ddH ₂ O, pH to 6.5 with concentrated NaOH |
| 100mM CaCl₂: | 5.55g CaCl ₂ in 500ml ddH ₂ O, |
| 1M Glycine: | 37.5g in 500ml ddH ₂ O, pH to 10.7 with concentrated NaOH |
| 10mM Oseltamivir Carboxylate: | 250mg GS4071 in 87.92ml ddH ₂ O, store at -80°C or 250mg Ro64-0802 in 64.7ml ddH ₂ O, store at -80°C |
| 10mM Zanamivir: | 200mg in 60.18ml ddH ₂ O, store at -80°C |
| 1mM MUNANA: | 25mg in 51ml MES assay buffer, store at -20C |

5.2 Working Solutions and Buffers

| | | |
|---------------------------------------|---|----------------|
| 100µM Oseltamivir Carboxylate: | 500µl of 10mM Oseltamivir carboxylate stock solution 49.5ml H ₂ O | Store at -20°C |
| 100µM Zanamivir: | 500µl 10mM Zanamivir stock solution 49.5ml H ₂ O | Store at -20°C |
| MES Assay Buffer: | 32.5mM MES: 50ml of 325mM MES stock solution 4mM CaCl ₂ : 20ml of 100mM CaCl ₂ stock solution ddH ₂ O: 430ml pH to 6.5 with concentrated NaOH | |
| 100µM MUNANA (per plate): | 300µl of 1mM stock solution 2.7ml MES assay buffer | |
| Stop Solution (500ml): | 0.1M Glycine: 50ml (1M stock) 25% Ethanol: 125ml absolute ethanol ddH ₂ O: 325ml pH to 10.7 with concentrated NaOH | |

6.0 NA Activity Determination (MUNANA Assay)

Use work instruction VW0574 for sample worksheet and results. Follow VW0771 for data management workflow.

The optimal virus sample dilution to standardise virus dose when measuring virus IC₅₀ to neuraminidase inhibitors (NIs) can be determined using this method.

Each assay should include subtype matched validated reference viruses. **Section 8** gives details of suitable references and validation criteria for the assay.

- 6.1 Add 20µl MES assay buffer to each well of a black 96 well flat bottomed plate.
- 6.2 Make duplicate two-fold dilutions of virus material, with a starting dilution of 1/2 by adding 20µl of the first virus to wells A1 and A2, 20µl of the second virus to wells A3 and A4 and so on until row A is filled (see table 4.1). Mix buffer and virus by pipetting up and down several times, taking care not to create aerosols.
- 6.3 Serial dilute the viruses down the plate by carrying over 20µl from row A to row B and so on, stopping at row G. Discard 20µl from row G. The final row of the plate contains buffer only as a blank control.
- 6.4 Prepare 3ml of MUNANA substrate working stock (100µM) per plate and add 30µl to each well including the blank row H.
- 6.5 Seal the plate(s) and incubate at 37°C for 60 minutes with shaking, in the dark.
- 6.6 Terminate the reaction by adding 150µl of stop solution to all wells.

Table 4.1 Plate layout for virus addition

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---------------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|
| 1/2 A | Virus 1 | Virus 1 | Virus 2 | Virus 2 | Virus 3 | Virus 3 | Virus 4 | Virus 4 | Virus 5 | Virus 5 | Virus 6 | Virus 6 |
| 1/4 B | | | | | | | | | | | | |
| 1/8 C | | | | | | | | | | | | |
| 1/16 D | | | | | | | | | | | | |
| 1/32 E | | | | | | | | | | | | |
| 1/64 F | | | | | | | | | | | | |
| 1/128 G | | | | | | | | | | | | |
| Buffer Only H | | | | | | | | | | | | |

- 6.7 Read the plate within 30 minutes of adding stop solution using the MUNANA test protocol on the Flurostar optima plate reader (appendix 2).
- 6.8 The data are plotted as relative fluorescence units (RFU) against virus dilution, with the mean blank (buffer only) value subtracted. This plot should yield a sigmoid dose-response curve (see section 8 and useful note 9.4).
- 6.9 To calculate standard virus dose for IC₅₀ testing, transfer the 96 well plate data to the MUNANA Results excel template (VW0574).
- 6.10 The standard virus dose is calculated by defining the virus dilution in which enzyme activity for a given isolate yields the equivalent level of fluorescence in one hour as 10µM of 4-methylumbelliferone sodium salt (useful note 9.1 and appendix 1). This cut off should be within the linear range of the enzyme activity curve. Examples of expected curves are given in appendix 3.

Appendix 1 gives details of how to generate the standard curve of 4-methylumbelliferone sodium salt.

7.0 Neuraminidase Inhibition Assay

Use work instruction VW0575 for sample worksheet and results. Follow VW0771 for data management workflow.

This section describes how to determine the IC_{50} of a virus to a neuraminidase inhibitor (NI).

For best results calculate the standard virus dose by measuring the NA activity of each virus (**NA Activity assay: section 6.0**) performing both assays on the same day (useful note 9.5).

Each assay should include the subtype matched validated reference viruses (**see section 8**).

- 7.1 Dilute each virus appropriately in MES assay buffer.
- 7.2 Add 10 μ l of diluted virus to 2 columns (wells A-G) of a black flat bottomed 96 well plate (i.e. column 1 +2 wells A-G virus 1, column 3+4 wells A-G virus 2 etc. See table 5.1).

Table 5.1 Plate layout for virus addition

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | |
|-------------|---|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|
| 1000 | A | Virus 1 | Virus 1 | Virus 2 | Virus 2 | Virus 3 | Virus 3 | Virus 4 | Virus 4 | Virus 5 | Virus 5 | Virus 6 | Virus 6 |
| 100 | B | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ |
| 10 | C | | | | | | | | | | | | |
| 1 | D | | | | | | | | | | | | |
| 0.1 | E | | | | | | | | | | | | |
| 0.01 | F | | | | | | | | | | | | |
| VC | G | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ | ↓ |
| Buffer Only | H | → | | | | | | | | | | | |

- 7.3 In an 8 well reservoir, prepare ten-fold dilutions of drug (see table 5.2)
- 7.4 Add 10µl of each drug dilution to a full row of a 96 well plate (i.e. Row A 1-12: 1,000nM, row B 1-12: 100nM, row C 1-12: 10nM etc). Ensure that the virus and drug are mixed.
- 7.5 Seal the plate(s) and incubate for 30 minutes at 37°C with shaking.
- 7.6 Prepare 3ml of MUNANA working stock (100µM) per plate and add 30µl of substrate to each well including the blank row H, ensuring virus/drug and substrate mix.
- 7.7 Seal the plate(s) and incubate at 37°C for 60 minutes with shaking, in the dark.
- 7.8 Terminate the reaction by adding 150µl of stop solution to all wells.

Table 5.2 Method of preparing drug dilutions*

| Step | Dilution series | Drug Concentration (nM) | 'In Assay' Concentration (nM) |
|------|---|--------------------------|-------------------------------|
| 1 | 100µl of 100µM working stock+1900µl MES | 5000 | 1000 |
| 2 | 150µl of step 1 +1350 MES | 500 | 100 |
| 3 | 150µl of step 2 +1350 MES | 50 | 10 |
| 4 | 150µl of step 3 +1350 MES | 5 | 1 |
| 5 | 150µl of step 4 +1350 MES | 0.5 | 0.1 |
| 6 | 150µl of step 5 +1350 MES | 0.05 | 0.01 |
| 7 | Buffer only | Virus/Substrate control | 0 |
| 8 | Buffer only | Substrate/Buffer control | 0 |

*NB: The drug dilution range can be changed to suit the known or suspected IC₅₀ values of the test isolates. A narrow range for more precise calculation of low IC₅₀ values (250nM-0.08nM) and an extended range for highly resistant isolates (10,000nM-0.1nM) have been calculated in work instruction VW0770.

7.9 The data are plotted as RFU against NA inhibitor concentration, with the mean blank (buffer only) value subtracted. This plot should yield a sigmoid dose-response curve (see section 8 and useful note 9.4).

7.10 To calculate the IC₅₀ values, transfer the plate data to the NAI Results excel template (VW0575).

7.11 IC₅₀ values are calculated for the duplicates independently, and the mean IC₅₀ taken as the final value. The RFU given by 50% of the virus control value is calculated, and the drug dilution corresponding to this level of fluorescence is the IC₅₀ value. An example is given in appendix 4.

7.12 Trends in IC₅₀ values for given subtypes and seasons are monitored, refer to SOP V6816 for details.

8.0 Neuraminidase Inhibition Assay Validation and Reference Criteria

- 8.1 Current reference strains and batches in use and the IC₅₀ and MUNANA validation limits are given in VW0772.
- 8.2 Reference viruses are included in all NA activity and IC₅₀ assays and should be subtype matched to the samples undergoing testing, wherever possible. If the sample subtype is unknown, reference strains of all subtypes should be used. Suggestions of reference suitable strains are given in the table below (Useful note 9.6).
- 8.3 All curves in both the NA activity and IC₅₀ assays should be manually checked for points which do not fit the sigmoidal shape and to ensure replicate curves match.
- 8.4 In the NA activity assay, the virus dose calculated from the replicates should be no more than one dilution factor apart.
- 8.5 The final IC₅₀ is a mean of independently calculated values from the replicates.
- 8.6 Validation limits for each reference virus should be determined. Maximum limits are defined as 3 standard deviations above and below the median IC₅₀. The median is calculated from a **minimum** of ten independent assays of the reference virus.
- 8.7 If a reference virus IC₅₀ value fails to meet validation criteria for a given drug, the test using that drug is invalidated and all samples repeated. Trends in reference IC₅₀ performance should be monitored. A batch which fails to meet validation criteria in three consecutive assays should be discarded.

9.0 Useful Notes


- 9.1 Fluorescence is measured in relative fluorescence units (RFU). Different fluorimeters have different ranges of values. Raw data values measured in RFU cannot be compared from one machine to another, and can only be compared from one assay to the next if settings are not changed.

- 9.2 Where biosafety levels dictate, plates can be read whilst sealed with an optical clear plate seal. This will reduce the overall fluorescence units slightly, but will not impact on the shape of the curves, nor the virus dilution/IC₅₀ value. In this instance, stop solution should be added to the plates, and then left to equilibrate to room temperature for 10 minutes prior to sealing the plate to minimise condensation.
- 9.3 Samples with low titre may have insufficient NA activity for inhibition testing and may give inaccurate IC₅₀ values. Only samples with HA titres of ≥ 16 and/or peak NA activity equivalent to 10 μ M 4-methylumbelliferone sodium salt can be reliably tested. Samples not meeting these criteria should be passaged to yield a higher titre.
- 9.4 Plotted data should yield sigmoid curves. Strains with low neuraminidase activity (despite reasonable HA titre) or flattened curves in presence of drug may be exhibiting reduced sensitivity to neuraminidase inhibitors, even if their IC₅₀ value is within the normal range. These samples should be subject to further characterisation.
- 9.5 The NA activity of influenza viruses, particularly those with mutations in the NA gene causing resistance can be unstable. Isolates should be stored at -80°C and kept at +4°C for minimal times. Virus dilutions calculated by MUNANA assay are only valid up to 24 hours after testing if isolates are stored at +4°C and should be assayed again if isolates are frozen/left longer.
- 9.6 As described in section 8, subtype matched NI sensitive and resistant viruses should be included as references in all assays. If resistant viruses are not available, subtype matched sensitive strains can be used as references provided the performance of such viruses in IC₅₀ assays is well characterised, (e.g. evaluated in 10-20 independent assays. This will allow a median value for the IC₅₀ of neuraminidase inhibitor susceptibility for that particular reference virus to be determined. Assay performance can then be validated according to the criteria described in section 8. Whilst this approach to standardisation will not absolutely guarantee the ability of the test to determine neuraminidase resistance it will provide a means to ensure day to day variation is minimised.

10.0 Appendices


APPENDIX 1: 4-METHYLUMBELLIFERONE SODIUM SALT STANDARD CURVE

1. Dilute 4-methylumbelliferone sodium salt (4-MUSS) in water to 1mM concentration.
2. Serial dilute the 4-MUSS in 1/2 steps, in stop solution (used for NA activity and IC₅₀ tests). The 4-MUSS must be titrated in stop solution to ensure that the fluorophore is fluorescing (requires high pH).
3. Pipette 200µl of each dilution of 4-MUSS onto the same plates which are used in IC₅₀ testing (black, flat bottomed).
4. Measure the fluorescence activity of the 4-MUSS titration series. The volume of 200µl must be measured as this is equal to the final volume which is measured in the NA activity and IC₅₀ assays.
5. An example curve for the 4-MUSS titration is given below.
6. Determine the RFU generated by 10µM 4-MUSS.
7. The number of RFU given by 10µM 4-MUSS can then be applied to curves generated by viral titrations to determine standard dose for IC₅₀ testing.
8. For example, based on the curve below, a cut off of 38000 RFU would be applied to all virus titrations. The total number of RFU will be different on different fluorimeters (useful note 9.1)



APPENDIX 2: PLATE READER TEST PROTOCOL

Note: Gain settings are specific to the machine




APPENDIX 3: ANALYSES OF MUNANA RESULTS

| Virus | Seasonal H1N1 Sensitive | | Seasonal H1N1 Oseltamivir Resistant | | Flu B Sensitive | | Flu B Oseltamivir Resistant | | H1N1 2009 Sensitive | | H1N1 Oseltamivir Resistant | |
|-------|-------------------------|-------|-------------------------------------|-------|-----------------|-------|-----------------------------|-------|---------------------|-------|----------------------------|-------|
| 2 | 47044 | 46833 | 43411 | 43619 | 47026 | 48649 | 47335 | 47848 | 46387 | 44472 | 43293 | 43427 |
| 4 | 45984 | 46279 | 36239 | 33912 | 50775 | 50145 | 46586 | 46583 | 45778 | 45969 | 38722 | 39179 |
| 8 | 41221 | 36954 | 24199 | 22160 | 40454 | 38032 | 34489 | 32862 | 46731 | 46563 | 29529 | 27309 |
| 16 | 27688 | 25737 | 14115 | 12554 | 26047 | 22445 | 20149 | 19835 | 44679 | 44788 | 16343 | 15715 |
| 32 | 17135 | 14903 | 7792 | 6798 | 15133 | 13690 | 11452 | 11638 | 35847 | 33868 | 9894 | 9517 |
| 64 | 8642 | 8234 | 3936 | 3155 | 7413 | 6593 | 6067 | 6128 | 22881 | 20878 | 5553 | 5051 |
| 128 | 4599 | 4310 | 2248 | 2025 | 3977 | 3507 | 3309 | 3445 | 14136 | 12257 | 3066 | 2948 |
| Blank | 458 | 446 | 464 | 454 | 467 | 444 | 446 | 460 | 455 | 444 | 461 | 457 |

| Average-Blank | Seasonal H1N1 Sensitive | | Seasonal H1N1 Oseltamivir Resistant | | Flu B Sensitive | | Flu B Oseltamivir Resistant | | H1N1 2009 Sensitive | | H1N1 Oseltamivir Resistant | |
|---------------|-------------------------|----------|-------------------------------------|----------|-----------------|----------|-----------------------------|----------|---------------------|----------|----------------------------|----------|
| 0.301029996 | 46589.33 | 46378.33 | 42956.33 | 43164.33 | 46571.33 | 48194.33 | 46880.33 | 47393.33 | 45932.33 | 44017.33 | 42838.33 | 42972.33 |
| 0.602059991 | 45529.33 | 45824.33 | 35784.33 | 33457.33 | 50320.33 | 49690.33 | 46131.33 | 46128.33 | 45323.33 | 45514.33 | 38267.33 | 38724.33 |
| 0.903089987 | 40766.33 | 36499.33 | 23744.33 | 21705.33 | 39999.33 | 37577.33 | 34034.33 | 32407.33 | 46276.33 | 46108.33 | 29074.33 | 26854.33 |
| 1.204119983 | 27233.33 | 25282.33 | 13660.33 | 12099.33 | 25592.33 | 21990.33 | 19694.33 | 19380.33 | 44224.33 | 44333.33 | 15888.33 | 15260.33 |
| 1.505149978 | 16680.33 | 14448.33 | 7337.33 | 6343.33 | 14678.33 | 13235.33 | 11183.33 | 10997.33 | 35392.33 | 33413.33 | 9439.33 | 9062.33 |
| 1.806179974 | 8187.33 | 7779.33 | 3481.33 | 2700.33 | 6958.33 | 6138.33 | 5612.33 | 5673.33 | 22426.33 | 20423.33 | 5098.33 | 4596.33 |
| 2.10720997 | 4144.33 | 3855.33 | 1793.33 | 1570.33 | 3522.33 | 3052.33 | 2854.33 | 2990.33 | 13681.33 | 11802.33 | 2611.33 | 2493.33 |
| 0.301029996 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 |
| 0.602059991 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 |
| 0.903089987 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 |
| 1.204119983 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 |
| 1.505149978 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 |
| 1.806179974 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 |
| 2.10720997 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 | 35000 |

| Find Conc | Seasonal H1N1 Sensitive | | Seasonal H1N1 Oseltamivir Resistant | | Flu B Sensitive | | Flu B Oseltamivir Resistant | | H1N1 2009 Sensitive | | H1N1 Oseltamivir Resistant | |
|----------------|-------------------------|-------------|-------------------------------------|-------------|-----------------|-------------|-----------------------------|-------------|---------------------|-------------|----------------------------|-------------|
| 0.301029996 | | | 0.554219367 | | | | | | | | | |
| 0.602059991 | | | 0.621670279 | | | | 0.879059679 | 0.846208533 | | | 0.709050673 | 0.696511217 |
| 0.903089987 | 1.031357134 | 0.943327511 | | | 1.007549576 | 0.952865739 | | | | | | |
| 1.204119983 | | | | | | | | | 1.461410577 | | | |
| 1.505149978 | | | | | | | | | 1.514258732 | | | |
| 1.806179974 | | | | | | | | | | | | |
| 2.10720997 | | | | | | | | | | | | |
| Concentration | 10.7487295 | 8.776624362 | 4.18475733 | 3.582773613 | 10.17535513 | 8.97151399 | 7.569369036 | 7.01792193 | 32.67824557 | 28.93413992 | 5.117415418 | 4.971772146 |
| Dilution | 5.4 | 4.4 | 2.1 | 1.8 | 5.1 | 4.5 | 3.8 | 3.5 | 16.3 | 14.5 | 2.6 | 2.5 |
| Final Dilution | 5 | | 2 | | 5 | | 4 | | 15 | | 3 | |




THIS COPY IS UNCONTROLLED UNLESS THE DOCUMENT DISPLAYS A CONTROLLED COPY STAMP

APPENDIX 4: ANALYSES OF NA INHIBITION ASSAY RESULTS

| Strain Name | Seasonal H1N1 Sensitive | | Seasonal H1N1 Oseltamivir Resistant | | Flu B Sensitive | | Flu B Oseltamivir Resistant | | H1N1 2009 Sensitive | | H1N1 Oseltamivir Resistant | |
|---------------|-------------------------|-------|-------------------------------------|-------|-----------------|-------|-----------------------------|-------|---------------------|-------|----------------------------|-------|
| 1000 | 515 | 547 | 13598 | 13856 | 1827 | 1589 | 7760 | 7315 | 508 | 528 | 8463 | 9095 |
| 100 | 1069 | 1075 | 33951 | 34144 | 9724 | 10196 | 25728 | 26179 | 962 | 959 | 27950 | 27532 |
| 10 | 6195 | 5454 | 38990 | 38796 | 26434 | 25404 | 35927 | 36181 | 4129 | 4010 | 34008 | 32705 |
| 1 | 22318 | 22259 | 40284 | 40521 | 40162 | 34296 | 39326 | 40378 | 20767 | 18191 | 35511 | 35257 |
| 0.1 | 37496 | 36824 | 39830 | 39691 | 36946 | 35288 | 40963 | 41339 | 34293 | 38282 | 36742 | 35774 |
| 0.01 | 39719 | 40702 | 40671 | 40216 | 41378 | 37116 | 40032 | 39582 | 36520 | 39726 | 36377 | 35019 |
| Virus Control | 38052 | 39977 | 39393 | 39357 | 39347 | 35693 | 39559 | 41017 | 37024 | 38823 | 37357 | 37883 |
| Blank | 483 | 497 | 487 | 482 | 492 | 428 | 473 | 492 | 470 | 472 | 672 | 555 |

| Average-Blank | Seasonal H1N1 Sensitive | | Seasonal H1N1 Oseltamivir Resistant | | Flu B Sensitive | | Flu B Oseltamivir Resistant | | H1N1 2009 Sensitive | | H1N1 Oseltamivir Resistant | |
|---------------|-------------------------|----------|-------------------------------------|----------|-----------------|----------|-----------------------------|----------|---------------------|----------|----------------------------|----------|
| 3 | 14.75 | 46.75 | 13097.75 | 13355.75 | 1326.75 | 1088.75 | 7259.75 | 6814.75 | 7.75 | 27.75 | 7962.75 | 8594.75 |
| 2 | 568.75 | 574.75 | 33450.75 | 33643.75 | 9223.75 | 9695.75 | 25227.75 | 25678.75 | 461.75 | 458.75 | 27449.75 | 27031.75 |
| 1 | 5694.75 | 4953.75 | 38489.75 | 38295.75 | 25933.75 | 24903.75 | 35426.75 | 35680.75 | 3628.75 | 3509.75 | 33507.75 | 32204.75 |
| 0 | 21817.75 | 21758.75 | 39783.75 | 40020.75 | 39661.75 | 33795.75 | 38825.75 | 39877.75 | 20266.75 | 17690.75 | 35010.75 | 34756.75 |
| -1 | 36995.75 | 36323.75 | 39329.75 | 39190.75 | 36445.75 | 34787.75 | 40462.75 | 39083.75 | 33792.75 | 37781.75 | 36241.75 | 35273.75 |
| -2 | 39218.75 | 40201.75 | 40170.75 | 39715.75 | 40877.75 | 36615.75 | 39531.75 | 39081.75 | 36019.75 | 39225.75 | 35876.75 | 34518.75 |
| Virus Control | 37551.75 | 39476.75 | 38892.75 | 38856.75 | 38846.75 | 35192.75 | 39058.75 | 40516.75 | 36523.75 | 38322.75 | 36856.75 | 37382.75 |
| 50% Cut | 18775.88 | 19738.38 | 19446.38 | 19428.38 | 19423.38 | 17596.38 | 19529.38 | 20258.38 | 18261.88 | 19161.38 | 18428.38 | 18691.38 |
| 3 | 18775.88 | 19738.38 | 19446.38 | 19428.38 | 19423.38 | 17596.38 | 19529.38 | 20258.38 | 18261.88 | 19161.38 | 18428.38 | 18691.38 |
| 2 | 18775.88 | 19738.38 | 19446.38 | 19428.38 | 19423.38 | 17596.38 | 19529.38 | 20258.38 | 18261.88 | 19161.38 | 18428.38 | 18691.38 |
| 1 | 18775.88 | 19738.38 | 19446.38 | 19428.38 | 19423.38 | 17596.38 | 19529.38 | 20258.38 | 18261.88 | 19161.38 | 18428.38 | 18691.38 |
| 0 | 18775.88 | 19738.38 | 19446.38 | 19428.38 | 19423.38 | 17596.38 | 19529.38 | 20258.38 | 18261.88 | 19161.38 | 18428.38 | 18691.38 |
| -1 | 18775.88 | 19738.38 | 19446.38 | 19428.38 | 19423.38 | 17596.38 | 19529.38 | 20258.38 | 18261.88 | 19161.38 | 18428.38 | 18691.38 |
| -2 | 18775.88 | 19738.38 | 19446.38 | 19428.38 | 19423.38 | 17596.38 | 19529.38 | 20258.38 | 18261.88 | 19161.38 | 18428.38 | 18691.38 |
| Virus Control | 18775.88 | 19738.38 | 19446.38 | 19428.38 | 19423.38 | 17596.38 | 19529.38 | 20258.38 | 18261.88 | 19161.38 | 18428.38 | 18691.38 |

| FINDIC50 | Seasonal H1N1 Sensitive | | Seasonal H1N1 Oseltamivir Resistant | | Flu B Sensitive | | Flu B Oseltamivir Resistant | | H1N1 2009 Sensitive | | H1N1 Oseltamivir Resistant | |
|-----------|-------------------------|-------------|-------------------------------------|-------------|-----------------|-------------|-----------------------------|-------------|---------------------|--------------|----------------------------|------------|
| 3 | | | 2.68807424 | 2.700678973 | | | 2.317140194 | 2.287339642 | | | 2.462943244 | 2.45237159 |
| 2 | | | | | 1.389609515 | 1.480495463 | | | | | | |
| 1 | 0.188666811 | 0.120224636 | | | | | | | 0.12049976 | | | |
| 0 | | | | | | | | | | -0.073198198 | | |
| -1 | | | | | | | | | | | | |
| -2 | | | | | | | | | | | | |
| RES | 0.188666811 | 0.120224636 | 2.68807424 | 2.700678973 | 1.389609515 | 1.480495463 | 2.317140194 | 2.287339642 | 0.12049976 | -0.073198198 | 2.462943244 | 2.45237159 |
| IC50 | 1.544069382 | 1.318938773 | 487.6118369 | 501.9713989 | 24.52502819 | 30.23398985 | 207.5583424 | 193.7936942 | 1.319774581 | 0.844893175 | 290.3643167 | 283.381562 |
| Mean IC50 | 1.4 | | 494.8 | | 27.4 | | 200.7 | | 1.1 | | 286.9 | |



THIS COPY IS UNCONTROLLED UNLESS THE DOCUMENT DISPLAYS A CONTROLLED COPY STAMP

11.0 SUMMARY OF REVISIONS

| | | | | |
|---------------------|-----|--------------------------|----|-------------------------------------|
| RETRAINING REQUIRED | YES | <input type="checkbox"/> | NO | <input checked="" type="checkbox"/> |
| N/A. First Issue. | | | | |

FOR INFORMATION ONLY